

PRODUCT CATALOG


Medicinal Maggots™ (maggot-impregnated gauze)

Medicinal Maggots™

Disinfected *Phaenicia* (=Lucilia) *sericata* larvae

Supplied as **maggot-impregnated gauze**, 2 individually wrapped pads per vial. Each vial contains 250 - 500 larvae (~200 larvae per gauze pad).

Catalog Number: MM250 \$295.00



Maggot MegaPak

Maggot MegaPak offers customized dressing options and everything you need, at one low single price. For simplicity, we have prepared the following complete treatment package for you, including shipping!


Items Included:

Medicinal Maggots™, 1 vial
Dressing of your choice
Ex: LeSoc™
Sheer Comfort™
LeGlue, if needed
Hydrocolloid Pad
Sterile gloves (1 pair)
Non-sterile gloves (2 pair)
Skin protectant (2)
Tape (1 roll)
Gauze, 1 pack

Includes FedEx Shipping, Standard Overnight

Catalog Number: MegaPk \$410.00

Maggot Dressings, à la carte



Sheer Comfort™

Polyester net dressing to keep the maggots on the wound. Using adhesive glue and tape, bond to a hydrocolloid pad or other foundation on the skin to build a maggot cage around the wound.

<u>Size:</u>	<u>Catalog Number:</u>	<u>Price:</u>
4 inches x 4 inches	CC44	\$10.00
8 inches x 8 inches	CC88	\$20.00
12 inches x 12 inches	CC1212	\$28.00
24 inches x 24 inches	CC2424	\$69.00
Custom sizes		Call for pricing
Nylon stocking, sterile knee-high	CCKH	\$16.00
Nylon pantyhose, sterile	CCPH	\$21.00

Maggot Dressings, à la carte



LeSoc™

Stocking-like dressing constructed from fixed-weave polyester monofilament netted fabric (like our Sheer Comfort™). Use to dress wounds with complicated 3-dimensional surfaces, like toes, feet, stumps. Similar to the Sterile Nylon Stockings, but the fixed-weave fabric better prevents against escaping maggots. Not for areas that may swell or become constricted.

<u>Size:</u>	<u>Catalog Number:</u>	<u>Price:</u>
1.75 inches x 2.5 inches LeSoc™ - Finger-sized Length: 2.5" / 6.4 cm; Width: 1.75" / 4.5 cm; Circumference: 3.5" / 8.9 cm; Diameter: 1" / 2.5 cm	LS1X2	\$7.00
3 inches x 4 inches LeSoc™ - Two-Finger sized Length: 4" / 10 cm; Width: 3" / 7.5 cm; Circumference: 6" / 15 cm; Diameter: 1.9" / 4.5 cm	LS3X4	\$12.00
4 inches x 6 inches LeSoc™ - Four-Finger-sized Length: 6" / 15 cm; Width: 4" / 10 cm; Circumference: 8" / 20 cm; Diameter: 2.5" / 6.4 cm	LS4X6	\$16.00
6 inches x 8 inches LeSoc™ - Hand-sized Length: 8" / 20 cm; Width: 6" / 15 cm; Circumference: 12" / 30 cm; Diameter: 3.8" / 9.5 cm	LS6X8	\$22.00
8 inches x 12 unches LeSoc™ - Anterior Foot Length: 12" / 30 cm; Width: 8" / 20 cm; Circumference: 16" / 40 cm; Diameter: 5" / 12.5 cm	LS8X12	\$31.00
12 inches x 18 inches LeSoc™ - Boot-sized Length: 18" / 45 cm; Width: 12" / 30 cm; Circumference: 24" / 60 cm; Diameter: 7.5" / 19 cm	LS12X18	\$46.00
18 inches x 24 inches LeSoc™ - Leg-sized Length: 24" / 60 cm; Width: 18" / 45 cm; Circumference: 36" / 91 cm; Diameter: 11.4" / 29 cm	LS18X24	\$62.00
24 inches x 36 inches LeSoc™ - Thigh-high Length: 36" / 91 cm; Width: 24" / 60 cm; Circumference: 48" / 122 cm; Diameter: 15" / 38 cm	LS24X36	\$84.00

Warning — since LeSoc™ is not elastic, care must be taken to ensure that it does not constrict or prevent adequate blood flow.

Maggot Dressings, à la carte



Hydrocolloid Pad

Used to make custom maggot therapy dressings at the bedside. Simply cut a hole in the pad to match the wound shape, place remaining pad as a foundation around the wound, place the maggots on the wound bed, and glue your porous cover (i.e., Sheer Comfort™ polyester net) to the hydrocolloid foundation.

Catalog Number: Price:

Hydrocolloid pad, 4 inches x 4 inches	HC44	\$14.00
Hydrocolloid pad, 6 inches x 6 inches	HC66	\$26.00



Le-Glu Liquid Adhesive By Nu-Hope

For gluing polyester net to hydrocolloid pad or other foundation. Shake (the bottle) well just before use.

Catalog Number: Price:

Liquid adhesive, 1	NHLG	\$16.00
--------------------	------	---------



Durapore Tape

Tape a window frame over the polyester net to secure it.

Catalog Number: Price:

1 inch wide x 10 yards; 1 roll	Tape	\$4.00
-----------------------------------	------	--------

Maggot Dressings, à la carte



Transparent membrane dressings

Used optionally to secure the maggot dressing to the skin, especially for wounds over areas prone to excessive movement (i.e., joints), perspiration or soiling (i.e., sacrum, coccyx). Placed as a frame over the net-hydrocolloid bond and extending onto the peripheral skin, can make maggot dressings more durable in areas predisposed to soiling and loosening (i.e., wounds over sacrum, coccyx, joints).

	<u>Catalog Number:</u>	<u>Price:</u>
~2 inches x 3 inches (brand may vary)	Tegaderm	\$1.50



Skin protective wipe

Used optionally to protect the healthy skin around the wound from rash, maceration or damage from the adhesives and wound drainage (assorted brands).

	<u>Catalog Number:</u>	<u>Price:</u>
Skin protective wipe (brand may vary)	SknPrep	\$0.76 each

Special Needs?

Do you need something special? Something you do not see in this catalog? Just ask.

We have more products than we have room in this catalog – from veterinary supplies to a variety of confinement dressings. We can even make custom orders.

Call: 949-679-3000.

Maggot Therapy for Animals?

Veterinary use of maggot therapy sometimes requires special modifications and materials.

Call us for veterinary products and services: 949-679-3000.



Disinfection processing under sterile hood

Monarch Labs' insectary--where it all begins!



At Monarch Labs, producing and delivering safe, effective and economical therapies for wound care patients is our highest priority.

We are the leading living medicine company in the country. We manufacture Medicinal Maggots™ (disinfected *Phaenicia sericata* larvae) in our production facility in Irvine, California. We distribute throughout North America via overnight courier.

Medicinal Maggots™ are regulated by the Food and Drug Administration as a *prescription only* medical device for treating pressure ulcers, venous stasis ulcers, neuropathic foot ulcers and non-healing traumatic or post surgical wounds.

Medicinal Maggots™ are highly perishable and require special dressings, which we also sell. They have been used for about 30 years to treat human and animal wounds.

Medicinal Maggots™ were the first living animals to receive marketing clearance by the FDA. They are reimbursed by major insurance companies. Since we began operations in 2005, over 2,000 health care institutions have ordered our products.

The Monarch team regards great customer service as fundamental to a successful business. Feel free to telephone or e-mail us with any questions or suggestions.

On behalf of the entire Monarch team,



James Kuo, MD, MBA
Manager / Monarch Labs
(949) 679-3000
jkuo@monarchlabs.com